

Republic of the Philippines  
**VISAYAS STATE UNIVERSITY**  
Visca, Baybay City, Leyte  
**SUPPLY & PROPERTY MANAGEMENT**

## INSPECTION & ACCEPTANCE REPORT

Supplier: **AG DEQUILLA TRADING**  
P.O No.: **PO-GF-CO-2025-06-0365**  
Office/Dept: **DBS**

IAR No: 25-07-0920  
Date: 2025-07-28  
66,699.90 GF-2023-08-01791

No.	UNIT	DESCRIPTION	Quantity		Date Received	Balance	Remarks
			Per P.O	Actual			
1	lot	2 Unit of Digital Electric Stereo Microscope with Camera Lucida	1	1	2025-07-21 #: DR-0135	0	Completed

**Specification:**

**Digital Electric Stereo Microscope with Camera Lucida**

**Hardware:**

- **Model:** AR tilting trinocular tube SZX2-ARTTR
- **Angle of Observation Tube:** 5 to 45 degrees
- **Maximum Power Consumption:** 10 W
- **Interpupillary Distance Adjustment Range:** 57 to 80 mm
- **Eyeiece Magnification:** 1.25X
- **Camera Magnification:** 1X
- **Eyeiece Clamping Knob:** Present
- **Light-Path Switching Mechanism:** Present
- **Front Panel Functions:** AR image brightness adjustment (7 levels), AR image ON/OFF
- **Input Connectors:** HDMI x 1, USB 2.0 (Type-C) x 1, DC jack x 1
- **Drive Voltage:** AC 100–240 V (AC adaptor)
- **Camera Lucida Adapter:** Present
- **Camera Lucida Included:** Yes

**Zoom Magnification Sensor:**

- **Main Functions:**
  - Set the zoom magnification (at click positions only)
  - Move slides created on the software forward and backward
- **Output Connector:** USB 2.0 Type-C

**Recording:**

- **File Saving Formats:** MP4, MOV
- **Video Codec:** H264
- **Audio Format:** MP3
- **Audio Codec:** MP3

**Additional Features:**

- **AR Image Recording:** The camera image can be saved with the AR image simultaneously
- **Recording Time:** Approximately one hour
- **Software Language:** English
- **Camera Control Functions:**
  - Exposure control: Switching between auto and manual
  - ISO Sensitivity Adjustment: ISO100, 200, 400
  - White Balance Adjustment Function: Equipped (one-touch adjustment)
- **Snapshot Acquisition:**
  - File Saving Formats: BMP, JPEG, PNG
  - Resolution:
    - DP23: 3088 × 2076
    - DP28: 4104 × 2174
  - The camera image can be saved with the AR image simultaneously

This microscope setup is designed for detailed imaging and recording with a high degree of flexibility and ease of use.

#### INSPECTION

Date Inspected: 2/28/15

☐ Inspected, verified and found in order as to quantity and specifications

☐ Others: \_\_\_\_\_

\_\_\_\_\_  
\_\_\_\_\_

**NOREVE JEAN M. AGAD**  
Inspector

#### ACCEPTANCE

Date Recieved: 2/28/15

☒ Complete

☐ Partial

**VIVIAN V. BALBARINO**  
Head, Supply & Property Management